

# Recombinant Trypsin 100% Animal Origin Free (AOF) Level 3

Storage: Liquid form: -15 °C to -20 °C; Powder form: 2 °C to 8 °C

#### Description

Trypsin (23.4 kDa) is a serine endoprotease that cleaves amide and ester bonds of amino acids Arginine (Arg) and Lysine (Lys).

Richcore's recombinant trypsin products in liquid and powder form are expressed in *Pichia pastoris* using an artificially synthesized gene as genetic construct. The products are manufactured in Richcore's GMP certified facility, are 100 % Animal Origin Free (AOF) level 3, stable, and free from viral contaminants and contaminating proteases such as Carboxypeptidase A and Chymotrypsin. Both products are greater than 98% pure.

#### Applications

The products have been developed for use in biopharmaceutical manufacturing processes requiring GMP grade, 100 % AOF level 3, trypsin, specifically in:

- 1. **Insulin Manufacturing**: To convert proinsulin to insulin by cleaving Arg amino acids at positions B31 and B32 and Lys and Arg amino acids at positions B64 and B65
- 2. **Vaccine Manufacturing**: To passage cells and accentuate viral propagation via proteolytic action

Additionally, both products can be used in:

- 1. Cell culturing for detaching adherent cells
- 2. Proteomics to cleave specific amino acids
- 3. As a process aid in food processing
- 4. In manufacturing cosmetics

#### Key Parameters and Data

Parameter	Specifications (Liquid)	Specifications (Powder)
Appearance	Clear, colorless to pale greenish	White to off-white,
	brown, odorless liquid	odorless powder
Protein by UV280 method	> 60 mg/ml	> 0.6 mg/mg powder
Specific Activity (Protein by	>3800 USP Units/mg	> 3800 Units/mg
UV 280 Method)		
Purity by GPC-HPLC	> 98%	> 98%
TVAC	< 10 cfu/ml	< 10 <sup>2</sup> cfu/g
ТҮМС	< 100 cfu/ml	< 10 <sup>4</sup> cfu/g
Pathogens	Absent	Absent

PDS/TRYP/01PS

#### Richcore Enzymes Pvt. Ltd

Plot No. 204 and 237, Bommasandra - Jigani Link Road, K.I.A.D.B Industrial Area, Bangalore, Karnataka – 560105, INDIA; Phone: +91- 8110 – 419777; Fax: +91 - 8110 – 416834; Email: info@richcoreindia.com



## **Cell Culturing:** Tested for cell dissociation on CHO-K1, Vero, MRC-5, and MDCK cells

#### Handling and Guidelines for Use

#### **Storage and Handling**

#### Liquid form

- 1. Store at -15 °C to -25 °C upon receipt
- 2. Material can sustain upto 5 freeze-thaw cycles. However, it is recommended to prepare one time use aliquots stored at -15  $^{\circ}$ C to -25  $^{\circ}$ C
- 3. Material should be thawed at room temperature

#### **Powder form**

1. Store at 2 °C to 8 °C upon receipt

### **Guidelines for Use**

#### Cell Culturing:

- 1. Add an appropriate amount from the aliquots/master vial to prepare one time use working solution in appropriate buffer at pH 7.2-7.8. For example, 1X PBS OR Cell Culture Media without Fetal Bovine Serum. It is recommended to use buffers free of Ca<sup>2+</sup> and Mg<sup>2+</sup>containing EDTA for optimal results
- 2. Mix the solution well for 1-2 minutes. Do not agitate, do not form bubbles
- 3. Check pH of the working stock. If necessary, adjust the pH to 7.2- 7.8 with 1M Tris base or 1N HCl
- 4. Filter the solution under aseptic conditions
- 5. Do not store over 3 weeks at 2 °C to 8 °C
- 6. Store at -15 °C to -25 °C for prolonged storage upto 2 months
- 7. Working solutions can sustain upto 5 freeze-thaw cycles without losing potency for cell dissociation. However, it is best practice to minimize the number of freeze-thaw cycles

# The products are not for human or animal consumption, or use as a drug. Please refer to the MSDS of the products for detailed safety profile

#### **Richcore Enzymes Pvt. Ltd**

Plot No. 204 and 237, Bommasandra - Jigani Link Road, K.I.A.D.B Industrial Area, Bangalore, Karnataka – 560105, INDIA; Phone: +91- 8110 – 419777; Fax: +91 - 8110 – 416834; Email: info@richcoreindia.com